

# CAPNV

## Rapid Intervention Capnograph



Small and  
portable



High degree of  
protection



Sound and light  
double alarm



## Overview

The CAPNV veterinary end-respiratory capnograph is designed to provide continuous monitoring of the end-tidal CO<sub>2</sub> concentration (EtCO<sub>2</sub>), capnography, and respiration rate (RR) of animals.



## Features

**Power supply:** 2XAAA battery/lithium battery; **Small but not simple, plug and play;** **Intelligent data transfer.**

**Power On/Off Button**

Turn off CapnoCOP by pressing the button "⏻" for 3s;

Press the "⏻" button to mute the alarm;

**Configuration**

Press the "↑" "↓" button to configure CapnoCOP.

**Zero calibration**

Press and hold "⏻" "↓" together to trigger the Zero Calibration process.

**Airway Adapter**

Large animals  
Small animals

**Parameter Display**


EtCO2  
AwRR  
Capnogram

**Status Display**

Bluetooth  
Battery status



- High standard Accuracy CO2 measurement by dual wave length
- Splash proof and rugged design for use in harsh conditions
- 4 buttons for convenient settings adjustments, alarm pause and power on/off
- Adjustable Sweep speed
- Selectable units mmHg, kPa, %
- Apnea detection and adjustable
- Configurable alarms for Respiration rate alarm, EtCO2 high and low, clogged or missing adapter
- Low battery indicator
- Visual and audible (buzzer) alarm
- Real-time Capnogram display
- Up to 4 hours continuous use
- Large animals and small animals airway adapters
- Ultra lightweight
- No annual routine calibration required
- Support english, italian, french, hungarian, polish (other language upon request)
- Bluetooth connectivity to export real time data



# Specifications

## Mechanical and Environment

|                          |  |
|--------------------------|--|
| Dimensions               | 60 x 47 x 43 mm (2.36 x 1.85 x 1.69 inches)  |
| Weight                   | ~80 g (complete device with batteries), < 350grs with transport pouch.                               |
| Mechanical robustness    | Meets the shock and vibration requirements for transport of SS-EN ISO 21647:2004 clause 21.102       |
| Temperature and Humidity | Operating: -5 to 50°C, 10 to 90% RH, non-condensing<br>Storage: -30 to 70°C, <90% RH, non-condensing |
| Atmospheric pressure     | 70 - 120 kPa   |

## EtCO<sub>2</sub> Measurement

|  |  |
|--|--|
| Transducer Type                        | Mainstream CO <sub>2</sub> Monitor   |
| Measuring method                       | Dual-wavelength, non-dispersive infrared spectroscopy  |
| Initialization Time                    | Capnogram, displayed within few seconds at an ambient temperature of 25°C, and full specifications within 1-2 minutes especially if temperature below 15°C |
| CO <sub>2</sub> Measurement Range      | 0 to 99 mmHg, 0 to 9.9%, 0 to 9.9 kPa  |
| CO <sub>2</sub> Accuracy               | ±(0.2% + 2% of reading)  |
| Total system response time             | < 1s   |
| Calibration                            | No routine user calibration required. An airway adapter zero is required when changing to a new airway adapter   |
| Recovery time after defibrillator test | Unaffected   |

## Respiration Rate Measurement

|                           |  |
|---------------------------|--|
| Respiratory rate Range    | 0~150 Breaths Per Minute (BPM)   |
| Respiration Rate Accuracy | ±1BPM  |
| Breath detect             | Adaptive threshold, minimum 1 % CO <sub>2</sub> change                 |
| Airway Adapters           | < 5 cc dead space ( large animals), < 1 cc dead space( small animals ) |

## Animals Machin Interface

|                     |  |
|---------------------|--|
| Display             | Active OLED color graphic display, 65 000 colors, 128 x 128 pixels<br>Multi-languages(English, French, Italian, Hungarian, polish and more to be supported)                                      |
| Front Panel Buttons | 4 buttons: power on/off, alarm pause, up/down selection  |
| Alarms              | No Adapter, Clogged Adapter, No Breath Detected, High AwRR, Low ETCO <sub>2</sub> , High ETCO <sub>2</sub> , Temperature out of range, Atmospheric pressure is over operating range, Low Battery |





## Power Specifications

|                   |  |
|-------------------|--|
| External Supply   | None   |
| Power Consumption | TBD  |
| Batteries         | Two AAA Cell batteries (2x1.5VDC): Alkaline IEC:LR03 or Energizer Ultimate Lithium L92 batteries. Use of other Lithium batteries may present a risk of fire or explosion |
| Battery life time | Duracell Plus Alkaline: ~4 hours<br>Energizer Ultimate Lithium L92: ~7 hours   |

## Regulatory Compliance

|                      |  |
|----------------------|--|
| Certificate          | CE   |
| Standards compliance | IEC 60601-1:2012<br>IEC 60601-1-8:2012<br>IEC 60601-1-12:2014<br>ISO 80601-2-55:2018 |

## Interfaaces

| Part Number | Designation                               |   | Included with device  |
|-------------|---|---|---|
| CAPNV       | CAPNV module                              | 1 |   |
| M401D-AL    | Disposable<br>Large animal airway adapter | 1 |  |
| M401D-N     | Small animal airway adapter               | 1 |  |
| RIC-P       | Pouch                                     | 1 |  |
| RIC-L       | Lanyard                                   | 1 |  |
| RIC-LI      | Lithium battery                           | 1 | (contain)   |

\* The data is subject to change without notice. Please refer to the manual for the contraindications and precautions