


Multi-parameter **Module MMP12**

It is used to monitor the patient's blood oxygen, blood pressure, electrocardiogram, respiration, body temperature, heart rate, and pulse



Dimension : 140 mm × 85 mm × 25 mm

12 leads
ECG

SPO2

NIBP

HR

RR

PR

2×TEMP

Features

- > With pulse oxygen and pulse rate monitoring function
- > With monitoring function of systolic blood pressure, diastolic blood pressure and mean pressure
- > With 12-lead ECG, 3/5-lead ECG, 1 breathing, 2 body temperature monitoring functions
- > Provide three patient modes: adult, child and newborn mode
- > The module thickness is only 25mm, exquisite and compact
- > Design of integrated blood pressure airway system
- > The working status of the real-time transmission module: hardware status, software status and sensor status, the upper computer can alarm in time according to the information
- > When the perfusion index is as low as 0.075%, the blood oxygen monitoring is accurate and reliable, meeting the application of surgery and ICU
- > Both blood oxygen and blood pressure adopt advanced algorithms, with anti-motion interference and weak signal measurement performance
- > Blood pressure measurement has three modes: manual, automatic and continuous mode
- > Blood pressure measurement has hardware and software dual overvoltage protection functions
- > Double timeout protection for blood pressure measurement (module timeout protection, provide timing trigger port of host computer)
- > The measurement results of the ECG measurement part include heart rate, body temperature, respiration and ST segment off set values of I, II, and V1 channels
- > ECG measurement has diagnosis, monitoring, HARDEST and surgery modes
- > ECG diagnostic analysis results include arrhythmia, conduction block, myocardial infarction, STT changes, ventricular hypertrophy, atrial enlargement, electrical axis deviation and other heart diseases

Specifications

ECG		RR	
Range	0.15mV-5.5mV	Range	0~120rpm
Accuracy	Undefined	Accuracy	15-120rpm: ± 2 rpm or $\pm 2\%$; Undefined(<15rpm)
Resolution	2.36uV/LSB	Base resistance	500-2000 Ω
lead type	3 Lead:I or II or III 5Lead:I,II,III,AVR,AVL,AVF,V1 12lead: I,II,III,AVR,AVL,AVF,V1,V2,V3,V4,V5,V6	Variable resistance	0.2-3.0 Ω

HR		SPO2	
Range	adult: 15~300bpm Child/Newborn: 15~350bpm	Range	0~100%
Accuracy	± 1 bpm	Accuracy	$\pm 2\%$ (70%~100%)Undefined(0~69%)
Resolution	1bpm	Resolution	1%

NIBP		PR	
Pressure Range	0-300mmHg	Range	25~300bpm
Pressure Accuracy	± 2 mmHg or $\pm 1\%$ (Whichever is greater)	Accuracy	± 3 bpm
Resolution	1mmHg	Resolution	1bpm
Systolic Range	Adult : 40~270mmHg Pediatric : 40~235mmHg Neonate : 40~130mmHg	PI	
Distolic Range	Adult : 10~210mmHg Pediatric:10~200mmHg Neonate:10~90mmHg	Range	0~20%
Mean Range	Adult : 20~230mmHg Pediatric : 20~225mmHg Neonate : 20~100mmHg	Accuracy	Undefined
Accuracy	The mean deviation< ± 5 mmHg The standard deviation<8mmHg	Resolution	0.001%

TEMP	
Range	0-50°C
Accuracy	$\pm 0.1^\circ\text{C}$
Resolution	0.1°C

Electical Specification

Power supply	DC.12V $\pm 5\%$
Power consumption	≤ 6 W
Communication	TTL,USART
Temperature	Operating 0°C~ 70°C Storage -40°C~ 70°C


Compliance

Standard	IEC 60601-2-27: 2011	IEC 80601-2-30: 2018	DIN EN ISO 81060-1: 2012
	BS EN 1060-3:1997 +A2: 2009	ISO 80601-2-61: 2017	AAMI EC57:2012/(R)2020

Interfaces

Power and Communication Interface

Power and communication interface socket pin definition diagram (Note: 4 and 6 pins are the same input, 3 and 5 pins are the same input)



Note 1: TXD is the data sent by the M103 module to the host.


Note 2: J1 socket pin pitch is 2.54mm

Power and communication interface socket J1 pins are defined as follows:

Pin No	Signal	Description
1	TXD	UART Sending data from module to host
2	RXD	UART Receiving data from host to the module
3	DGND	12V battery ground
4	+12V	Power Supply input 12V DC
5	DGND	12V battery ground
6	+12V	Power Supply input 12V DC

Body temperature probe interface


The temperature probe interface J3 is defined as follows, and the socket spacing is 2.54mm:



Pin No	Signal	Description
1	TEMP1	Temperature sensor 1 positive direction input
2	TGND	Temperature sensor 1 negative input
3	TEMP2	Temperature sensor 2 positive direction input
4	TGND	Temperature sensor 2 negative input


ECG lead input interface

The 3/5-lead ECG lead input interface J4 is defined as follows, and the socket pin spacing is 2.54mm:




Pin No	Signal	Description
1	RL	RL(right leg)
2	ECG SHIELD	Signal shielding wire
3	V1	The chest electrode:4th Intercostal(IC)space at right border of sternum
4	LL	LL(left leg)
5	LA	LA(left arm)
6	RA	RA(right arm)

The 12-lead ECG lead input interface J4/J10 is defined as follows, and the socket pin spacing is 2.54mm:




A schematic diagram of the location of the lead electrodes is shown below:



Pin No	Signal	Description
1	RL	RL(right leg)
2	ECG SHIELD	Signal shielding wire
3	V1	The chest electrode:4th Intercostal(IC) Space at right border of sternum
4	LL	LL(left leg)
5	LA	LA(left arm)
6	RA	RA(right arm)
7	--	Reserved
8	V6	Left midaxillary line at the horizontal level of V4
9	V5	Left midaxillary line at the horizontal level of V4
10	V4	5th IC space on left midclavicular line
11	V3	Midway between V2 and V4
12	V2	The fourth intercostal space at the left edge of sternum

SpO2 probe interface














The pin definition of sensor interface socket J7 is as follows, and the socket spacing is 2mm:







Pin No	Signal	Description
1	1_WIRE	Probe identification signal line
2	GND	power ground
3	IR	IR LED positive electrode
4	RED	RED LED positive electrode
5	SPO2-	Light Signal Negative Electrode
6	SPO2+	Light Signal Positive Electrode

Purchase Guider

Code	Name	Quantity	Description	Image
042-130401-00	MMP12 Module	x 1		
024-030303-03	Internal tube for NIBP	x 1	Connect with M103 Module and Tee joint	
023-050103-00	Tee joint	x 1	Connect with Internal tube	
022-990101-00	Connector board	x 1	Connect with M103 Module	
041-004003-00	ECG Cable 10-lead defibrillation proof, European standard, button type	Optional	Connect with ECG interface	
041-004003-04	ECG Cable 10-lead defibrillation proof, American Standard, button type	Optional	Connect with ECG interface	
041-009004-00	SPO2 cable 6-pin SpO2 extension cable,DB9 type	x1	Connect with connector board SPO2 interface	
041-002004-01	SPO2 probe finger adult Applicable to adult	Optional	Connect with SPO2 cable	
041-002005-01	SPO2 probe finger child Applicable to Child	Optional	Connect with SPO2 cable	
041-002009-00	SPO2 probe wrap Applicable to adult /Pediatric/Infant	Optional	Connect with SPO2 cable	
041-002027-01	SPO2 soft probe finger adult Applicable to adult	Optional	Connect with SPO2 cable	
041-002028-01	SPO2 soft probe finger child Applicable to Child	Optional	Connect with SPO2 cable	
041-002030-00	Disposable SPO2 cable Non-adhesive	Optional	Connect with SPO2 cable	
041-002006-01	Disposable SPO2 cable Adhesive foam	Optional	Connect with SPO2 cable	
041-002005-01	Disposable SPO2 cable Adhesion	Optional	Connect with SPO2 cable	

Code	Name	Quantity	Description	Image
041-007002-00	NIBP External tube Reusable NIBP hose	x1	Connect with pressure sensor and cuff	
041-001024-01	Thigh Cuff Arm Range : 46-66cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001025-01	Large Adult Cuff Arm Range : 33-47cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001026-01	Adult Cuff Arm Range : 25-35cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001027-01	Adult Cuff Arm Range : 20-28cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001028-01	Child Cuff Arm Range : 18-26cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001029-01	Infant Cuff Arm Range : 10-19cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001032-01	Neonate Cuff Arm Range : 6-11cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001050-01	Neonate 1 size (disposable) Arm Range : 3-6cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001049-01	Neonate 2 size (disposable) Arm Range : 4-8cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001048-01	Neonate 3 size (disposable) Arm Range : 6-11cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001047-01	Neonate 4 size (disposable) Arm Range : 7-13cm	Optional	Connect with pressure sensor or NIBP external tube	
041-001046-01	Neonate 5 size (disposable) Arm Range : 8-15cm	Optional	Connect with pressure sensor or NIBP external tube	
041-003005-01	Temperature probe Adult Skin	Optional x 2	Connect with TEMP interface	

Code	Name	Quantity	Description	Image
041-003003-00	Temperature probe Pediatric General	Optional x 2	Connect with TEMP interface	
041-003007-00	Temperature probe Pediatric Skin(Disposable)	Optional x 2	Connect with TEMP interface	
041-003006-00	Temperature probe Adult General(Disposable)	Optional x 2	Connect with TEMP interface	
041-003008-00	Temperature probe Adult Skin(Disposable)	Optional x 2	Connect with TEMP interface	

* The data is subject to change without notice. Please refer to the manual for the contraindications and precautions