





Pulse Arterial Oxygen Saturation **Module MISPO2**

The MISPO2 module is used to monitor the patient's pulse oximetry, pulse, and perfusion index



Dimension : 88.9 mm × 50.8 mm × 7 mm



Features

- > With the monitoring function of pulse oxygen, pulse rate, perfusion index and variation index
- > Patient model: adult, child and newborn
- > Real-time transmission module working status: hardware status, software status and sensor status
- > It has the function of setting the average time of calculation parameters, and obtains the response time of different calculation parameters
- > When the blood perfusion index is as low as 0.075%, the measurement is accurate and reliable, which can meet the application of surgery and ICU
- > Adopt advanced anti-motion algorithm, with strong anti-motion performance

Specifications

SPO2

Range: 0~100%
 Accuracy: $\pm 2\%$ (70%~100%)
 Undefined (0~69%)
 Resolution: 1%

PR

Range: 25~300bpm
 Accuracy: ± 3 bpm
 Resolution: 1bpm

PI

Range: 0~20%
 Accuracy: Undefined
 Resolution: 0.001%

PVI

Range: 100%
 Accuracy: Undefined
 Resolution: 1%

Electical Specification

Power supply DC.12V \pm 10%

Power consumption ≤ 2 W

Communication TTL,USART

Temperature Operating 0°C~ 70°C Storage -40°C~ 70°C

Compliance

Standard ISO 80601-2-61: 2017
 IEC 60601-1-2:2014

Dimensions




Figure 1 Top view of SpO2 circuit board





Figure 2 Front view of SpO2 board

Interfaces

Power and Communication Interface

The pin definitions of the power and communication interface socket (J3) are as follows:




Note 1: TXD is the data interface sent by the LowPulseStr blood oxygen module to the host, in TTL level mode.

Note 2: The pin pitch is 2.0mm

Pin No	Signal	Description
1	TXD	UART Sending data from module to host
2	RXD	UART Receiving data from host to the module
3	GND	12V battery ground
4	+12V	Power Supply input 12V DC

SpO2 probe interface




Note 1: The pin pitch is 2.0mm

The pin definitions of the SpO2 interface socket (J2) are as follows:

Pin No	Signal	Description
1	1_WIRE	Probe identification signal line
2	GND	Power ground
3	IR	IR LED positive
4	RED	RED LED positive
5	SPO2-	Light Signal Negative Electrode
6	SPO2+	Light Signal Positive Electrode

Temperature probe interface



Note 1: Pin pitch is 2.54mm

The temperature probe interface socket (J4) pins are defined as follows:

Pin No	Signal	Description
1	TEMP1	Temperature sensor 1, positive terminal
2	TGND	Temperature sensor 1, reverse side
3	TEMP2	Temperature sensor 2, positive terminal
4	TGND	Temperature sensor 2, reverse side

Purchase Guider - 1

Code	Name	Quantity	Description	Image
022-040201-00	MISPO2 board	x1		
009-030108-00	Internal Cable with ODU connector	x1	Connect to M201B	
041-009001-00	External Cable ODU	Optional	Connect to ODU connector	
041-002004-01	SPO2 probe Finger adult Applicable to adult	Optional	Connect to External Cable	
041-002005-01	SPO2 probe Finger child Applicable to child	Optional	Connect to External Cable	
041-002009-00	SPO2 probe warp Applicable to Infant	Optional	Connect to External Cable	
041-002030-00	Disposable SPO2 cable Non-adhesive	Optional	Connect with SPO2 cable	
041-002006-01	Disposable SPO2 cable Adhesive foam	Optional	Connect with SPO2 cable	
041-007002-00	Disposable SPO2 cable Adhesion	Optional	Connect with SPO2 cable	

* The data is subject to change without notice. Please refer to the manual for the contraindications and precautions